


Date 8.1.09Course 2230-LName DR. STASSI DIMAGGIOSection 4Experiment MELTING POINT DETERMINATION


PURPOSE: TO DETERMINE THE IDENTITY OF AN UNKNOWN COMPOUND BASED ON ITS MELTING POINT

EQUATIONS: NONE


STRUCTURES:



NAPHTHALENE




ACETANILIDE




UREA




4-NITROTOLUENE




1-NAPHTHOL




2-METHYLBENZOIC ACID




BENZOIC ACID




SALICYLIC ACID



SULFANILAMIDE



SUCCINIC ACID



3,5-DINITROBENZOIC ACID

Date 8.1.09Course 2230-LName DR. STASSI DIMAGGIOSection 4Experiment MELTING POINT DETERMINATION

PROCEDURE: OBTAIN A SMALL AMOUNT OF SAMPLE IN A CAPILLARY TUBE (< 2mm HIGH) AND PACK TO THE BOTTOM OF THE TUBE. REPEAT FOR EACH KNOWN. PLACE SAMPLES IN MELTEMP WITH THERMOMETER AND BEGIN HEATING THE MELTEMP AT A SETTING WHERE THE TEMPERATURE RISES $\approx 1^{\circ}\text{C}/\text{MINUTE}$.

RECORD MELTING TEMPERATURE RANGE WHEN THE FIRST DROPS FORM TO WHEN THE SAMPLE IS LIQUID. COOL MELTEMP AND REPEAT FOR UNKNOWN. PERFORM A MIXED MELTING POINT FOR THE 3 MOST LIKELY CHOICES FOR THE UNKNOWN.

REFERENCE: CRAIG, K. et al. ORGANIC CHEMISTRY LABORATORY MANUAL, 6TH EDITION WILEY, NEW YORK, 2005, p 32,
WWW.XULA.EDU/CHEMISTRY/DEPARTMENT/ORGANIC/ORG0.HTML

DATA: ① A SAMPLE OF EACH OF NAPHTHALENE, ACETANILIDE AND UREA WERE OBTAINED
② THE SAMPLES WERE GROUND INTO A FINE POWDER

Date 8.1.09Course 2230LName DR. STASSI DIMAGGIOSection 4Experiment MELTING POINT DETERMINATION

- ③ ^{THE} SAMPLES WERE EACH COLLECTED IN A CAPILLARY TUBE BY TAPPING THE OPEN END OF THE TUBE ON THE COMPOUND AND TURNING THE TUBE OVER AND PACKING THE SAMPLE TO THE BOTTOM
- ④ A THERMOMETER WAS PLACED IN THE MELTEMP ALONG WITH THE 3 KNOWN SAMPLES AND THE DIAL WAS SET TO 3. (8:43am)
- ⑤ AFTER 5 MINUTES, THE TEMPERATURE STOPPED INCREASING AND THE DIAL WAS SET TO 4. (8:48am)
- ⑥ AT 8:52am THE DIAL WAS INCREASED TO 6 AND THE TEMP WAS INCREASING AT $\approx 1^\circ\text{C}/\text{MIN}$.

	OBSERVATION	1 ST DROP FORMED ($T^\circ\text{C}$)	COMPLETELY MELTED ($T^\circ\text{C}$)
ACETANILIDE	WHITE CRYSTALS	109.1 $^\circ\text{C}$	112.0 $^\circ\text{C}$
NAPHTHALINE	FLAT, WHITE FLAKES	79.5 $^\circ\text{C}$	81.5 $^\circ\text{C}$
UREA	OFF-WHITE POWDER	128.8 $^\circ\text{C}$	131.0 $^\circ\text{C}$
UNKNOWN # 342		98.0 $^\circ\text{C}$	101.5 $^\circ\text{C}$
UNK + 1-NAPHTHOL		82.5 $^\circ\text{C}$	87.6 $^\circ\text{C}$
UNK + 2-METHYL BENZOIC ACID		99.8 $^\circ\text{C}$	102.5 $^\circ\text{C}$
UNK + ACETANILIDE		95.5 $^\circ\text{C}$	100.0 $^\circ\text{C}$

Date 8.1.09Course 2230LName DR. STASSI DiMAGGIOSection 4Experiment MELTING POINT DETERMINATION

⑦ ONCE THE KNOWNs WERE MELTED AND THE TEMPERATURE RECORDED (SEE TABLE)

THE MELTEMP WAS COOLED TO 40.0°C

⑧ UNKNOWN # 342 WAS PREPARED ACCORDING TO PROCEDURE (STEPS 1-3) AND THE MELTING POINT RANGE WAS OBTAINED AND RECORDED

$98.0^{\circ}\text{C} - \cancel{101.0^{\circ}\text{C}} \quad 101.5^{\circ}\text{C}$
SD

⑨ SAMPLES OF PURE 1-NAPHTHOL, 2-METHYL BENZOIC ACID, AND ACETANILIDE WERE OBTAINED. SMALL AMOUNTS OF UNK # 342 WERE MIXED WITH EACH KNOWN AND THE MIXED MELTING POINT OF EACH WAS OBTAINED AND RECORDED (SEE DATA TABLE)

Stassi DiMaggio
8/1/09

